



Welltech Instrument Company Ltd

Temperature Chamber YDT Series

YDT Series



Welltech provide high-quality Temperature Chamber with different capacity that meet every customer's need, custom design also available.

Key Features:

- Temperature Range: -70 degC to 180degC
- LCD programmable controller, can program up to 120 pattern, 1000 segment
- Stability: +/-0.5degC
- Uniformity: +/-1 degC
- Time rise from -70 degC to 100degC is 30-45mins
- Time drop from 20 degC to -70degC is 30-45mins
- Material: All S.S.
- With RS232 to computer



Specification:

Model	Capacity	Temperature Range	Accuracy	Internal Dimension (HxWxD)cm
YDT63	63L	0 to 100 degC	+/-1 degC	45x40x35
YDT63J	63L	-20 to 100 degC	+/-1 degC	45x40x35
YDT63R	63L	-40 to 100 degC	+/-1 degC	45x40x35
YDT63Z	63L	-70 to 100 degC	+/-1 degC	45x40x35
YDT80	80L	0 to 100 degC	+/-1 degC	50x40x40
YDT80J	80L	-20 to 100 degC	+/-1 degC	50x40x40
YDT80R	80L	-40 to 100 degC	+/-1 degC	50x40x40
YDT80Z	80L	-70 to 100 degC	+/-1 degC	50x40x40
YDT120	120L	0 to 100 degC	+/-1 degC	60x50x40
YDT120J	120L	-20 to 100 degC	+/-1 degC	60x50x40
YDT120R	120L	-40 to 100 degC	+/-1 degC	60x50x40
YDT120Z	120L	-70 to 100 degC	+/-1 degC	60x50x40
YDT150	150L	0 to 100 degC	+/-1 degC	60x50x50
YDT150J	150L	-20 to 100 degC	+/-1 degC	60x50x50
YDT150R	150L	-40 to 100 degC	+/-1 degC	60x50x50
YDT150Z	150L	-70 to 100 degC	+/-1 degC	60x50x50
YDT225	225L	0 to 100 degC	+/-1 degC	75x50x60
YDT225J	225L	-20 to 100 degC	+/-1 degC	75x50x60
YDT225R	225L	-40 to 100 degC	+/-1 degC	75x50x60
YDT225Z	225L	-70 to 100 degC	+/-1 degC	75x50x60
YDT408	408L	0 to 100 degC	+/-1 degC	85x60x80
YDT408J	408L	-20 to 100 degC	+/-1 degC	85x60x80
YDT408R	408L	-40 to 100 degC	+/-1 degC	85x60x80
YDT408Z	408L	-70 to 100 degC	+/-1 degC	85x60x80
YDT800	800L	0 to 100 degC	+/-1 degC	100x100x80
YDT800J	800L	-20 to 100 degC	+/-1 degC	100x100x80
YDT800R	800L	-40 to 100 degC	+/-1 degC	100x100x80
YDT800Z	800L	-70 to 100 degC	+/-1 degC	100x100x80
YDT960	960L	0 to 100 degC	+/-1 degC	100x100x96
YDT960J	960L	-20 to 100 degC	+/-1 degC	100x100x96
YDT960R	960L	-40 to 100 degC	+/-1 degC	100x100x96
YDT960Z	960L	-70 to 100 degC	+/-1 degC	100x100x96
YDT1000	1000L	0 to 100 degC	+/-1 degC	100x100x100
YDT1000J	1000L	-20 to 100 degC	+/-1 degC	100x100x100
YDT1000R	1000L	-40 to 100 degC	+/-1 degC	100x100x100
YDT1000Z	1000L	-70 to 100 degC	+/-1 degC	100x100x100